

# NeuroCare 40Hz

## Personalized Neural Entrainment Therapy

### What It Does

Gently restores 40 Hz brainwave patterns disrupted in Alzheimer's and MCI patients using non-invasive sound stimulation — automatically adjusted to each patient's brain state.

### How It Works


1. Patient wears lightweight EEG headset during therapy session
2. System continuously reads brain activity in real time
3. AI predicts optimal stimulation timing 5–10 seconds ahead
4. 40 Hz audio stimulus activates when brain synchrony is lowest

### Why It Matters

- Non-invasive: no drugs, no implants, no side effects
- Personalized: each session adapts to the patient's unique response
- Evidence-based: validated on EEG recordings from 35 participants
- Proactive: predicts the right moment before the window passes

### Key Result

**72% targeting accuracy** vs 64% reactive baseline — statistically significant improvement across 35 subjects (Cohen's  $d = 1.31$ ,  $p < 0.001$ )



Scan to try the live demo



Share your feedback